

Product Datasheet

CHST1 peptide (orb218158)

Description	The peptide is used to block anti-CHST1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CHST1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214802.
Purity	> 85%
Source	Synthetic
Uniprot ID	O43916, Q5RJQ0, Q9EQC0
Entrez	76969, 8534, 295934
Expiration Date	12 months from date of receipt.