

Product Datasheet

p63 peptide (orb218165)

Description	p63 peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	TP63
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214809.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H3D4
Entrez	8626
Expiration Date	12 months from date of receipt.