

Product Datasheet

KAT2B peptide (orb218195)

Description	KAT2B peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	KAT2B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214839.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92831, Q9JHD1
Entrez	18519, 8850
Expiration Date	12 months from date of receipt.