

Product Datasheet

STAM2 peptide (orb218202)

Description	The peptide is used to block anti-STAM2
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	STAM2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214846.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5XHY7, O88811, O75886
Entrez	311030, 56324, 10254
Expiration Date	12 months from date of receipt.