

Product Datasheet

SOCS3 peptide (orb218205)

Description	The peptide is used to block anti-SOCS3 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	SOCS3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214849.
Purity	> 85%
Source	Synthetic
Uniprot ID	035718, 014543, 088583
Entrez	12702, 9021, 89829
Expiration Date	12 months from date of receipt.