

Product Datasheet

CCRL2 peptide (orb218206)

Description	The peptide is used to block anti-CCRL2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CCRL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214850.
Purity	> 85%
Source	Synthetic
Uniprot ID	000421, 035457
Entrez	54199, 9034
Expiration Date	12 months from date of receipt.