

Product Datasheet

PIAS2 peptide (orb218208)

Description	PIAS2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PIAS2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214852.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6AZ28, Q8C5D8, O75928
Entrez	83422, 17344, 9063
Expiration Date	12 months from date of receipt.