

Product Datasheet

TBX15 peptide (orb218214)

Description	The peptide is used to block anti-TBX18 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	TBX15; TBX18
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214858.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9EPZ6, O95935, Q96SF7, O70306
Entrez	9096, 76365, 21384
Expiration Date	12 months from date of receipt.