

Product Datasheet

AIFM1 peptide (orb218219)

Description	The peptide is used to block anti-AIFM1 antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	AIFM1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214863.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Z0X1 , Q9JM53 , O95831
Entrez	83533 , 26926 , 9131
Expiration Date	12 months from date of receipt.