

Product Datasheet

RPS6KA5 peptide (orb218239)

Description	The peptide is used to block anti-MSK1
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	RPS6KA5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214883.
Purity	> 85%
Source	Synthetic
Uniprot ID	075582, Q8C050
Entrez	9252, 73086
Expiration Date	12 months from date of receipt.