

Product Datasheet

UBE4A peptide (orb218247)

Description	The peptide is used to block anti-UBE4A antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	UBE4A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214891.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q14139 , Q6P7A2
Entrez	315608 , 9354
Expiration Date	12 months from date of receipt.