

Product Datasheet

EIF2AK3 peptide (orb218257)

Description	The peptide is used to block anti-PERK antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	EIF2AK3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214901.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NZJ5
Entrez	9451
Expiration Date	12 months from date of receipt.