

Product Datasheet

TOMM20 peptide (orb218289)

Description	The peptide is used to block anti-TOM20 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TOMM20
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214933.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9DCC8 , Q62760 , Q15388
Entrez	67952 , 266601 , 9804
Expiration Date	12 months from date of receipt.