

## Product Datasheet

### ZC3H11A peptide (orb218296)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-ZC3H11A antibody reactivity. |
| <b>Reactivity</b>          | Human, Mouse, Rat                                              |
| <b>Tested Applications</b> | Blocking                                                       |
| <b>Target</b>              | ZC3H11A                                                        |
| <b>Form/Appearance</b>     | Lyophilized powder                                             |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                   |
| <b>Note</b>                | For research use only                                          |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb214940.       |
| <b>Purity</b>              | > 85%                                                          |
| <b>Source</b>              | Synthetic                                                      |
| <b>Uniprot ID</b>          | <b>O75152, Q6NZF1</b>                                          |
| <b>Entrez</b>              | <b>70579, 9877</b>                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                |