

Product Datasheet

OR5M3 Antibody - middle region (orb2182965)

Description	OR5M3 Antibody - middle region
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthesized peptide derived from the Internal region of human Olfactory receptor 5M3. Location: 30-110aa.
Preservatives	Liquid. 1 mg/mL in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid. 1 mg/mL in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Storage	Store at -2°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	Application Info: WB 1:500-1:2000; IF 1:200-1:1000; ELISA 1:5000;
Isotype	IgG
Clonality	Polyclonal
MW	35 kDa
Uniprot ID	Q8NGP4
Expiration Date	12 months from date of receipt.