

Product Datasheet

G3BP2 peptide (orb218300)

Description	The peptide is used to block anti-G3BP2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	G3BP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214944.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UN86, P97379
Entrez	23881, 9908
Expiration Date	12 months from date of receipt.