

Product Datasheet

FGF19 peptide (orb218304)

Description	The peptide is used to block anti-FGF19 antibody reactivity.
Reactivity	Human, Mouse, Primate
Tested Applications	Blocking
Target	FGF19
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214948.
Purity	> 85%
Source	Synthetic
Uniprot ID	095750
Entrez	9965
Expiration Date	12 months from date of receipt.