

Product Datasheet

HDAC5 peptide (orb218309)

Description	The peptide is used to block anti-Histone Deacetylase 5 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	HDAC5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214953.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Z2V6 , Q9UQL6
Entrez	10014
Expiration Date	12 months from date of receipt.