

Product Datasheet

BCKDK peptide (orb218330)

Description	The peptide is used to block anti-BCKDK antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	BCKDK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214974.
Purity	> 85%
Source	Synthetic
Uniprot ID	O14874, O55028, Q00972
Entrez	29603, 10295, 12041
Expiration Date	12 months from date of receipt.