

Product Datasheet

SIRPB1 peptide (orb218332)

Description	The peptide is used to block anti-CD172b antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	SIRPB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214976.
Purity	> 85%
Source	Synthetic
Uniprot ID	000241
Entrez	10326
Expiration Date	12 months from date of receipt.