

Product Datasheet

AKR1A1 peptide (orb218333)

Description	The peptide is used to block anti-AKR1A1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	AKR1A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214977.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JII6, P51635, P14550
Entrez	58810, 10327, 78959
Expiration Date	12 months from date of receipt.