

Product Datasheet

PIAS3 peptide (orb218338)

Description	The peptide is used to block anti-PIAS3 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	PIAS3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214982.
Purity	> 85%
Source	Synthetic
Uniprot ID	O70260 , O54714 , Q9Y6X2
Entrez	10401 , 83614 , 229615
Expiration Date	12 months from date of receipt.