

Product Datasheet

YAP1 peptide (orb218342)

Description	The peptide is used to block anti-YAP1 antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	YAP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214986.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q2EJA0 , P46937 , P46938
Entrez	22601 , 363014 , 10413
Expiration Date	12 months from date of receipt.