

Product Datasheet

PGRMC2 peptide (orb218343)

Description	The peptide is used to block anti-PGRMC2 antibody reactivity.
Reactivity	Human, Porcine, Rat
Tested Applications	Blocking
Target	PGRMC2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214987.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5XIU9 , O15173
Entrez	10424 , 361940
Expiration Date	12 months from date of receipt.