

Product Datasheet

MRSA peptide (orb218351)

Description	MRSA peptide
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	SLC9A6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214995.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92581
Entrez	10479
Expiration Date	12 months from date of receipt.