

Product Datasheet

KSR1 Antibody (Phospho-Ser392) (orb2183542)

Description	KSR1 Antibody (Phospho-Ser392)
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthesized peptide derived from human Ksr-1 around the phosphorylation site of S392. Location: 330-410aa.
Preservatives	Liquid. 1 mg/mL in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid. 1 mg/mL in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Storage	Store at -2°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	Application Info: WB 1:500-1:2000; IHC 1:100-1:300; ELISA 1:20000;
Isotype	IgG
Clonality	Polyclonal
MW	115 kDa
Uniprot ID	Q61097, Q8IVT5
Expiration Date	12 months from date of receipt.