

Product Datasheet

DCTN2 peptide (orb218355)

Description	The peptide is used to block anti-Dynactin p50 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	DCTN2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214999.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13561 , Q6AYH5 , Q99KJ8
Entrez	69654 , 299850 , 10540
Expiration Date	12 months from date of receipt.