

Product Datasheet

GPR75 peptide (orb218384)

Description	The peptide is used to block anti-GPR75 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GPR75
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215028.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6X632, O95800
Entrez	10936, 237716
Expiration Date	12 months from date of receipt.