

Product Datasheet

CDC37 peptide (orb218396)

Description	The peptide is used to block anti-CDC37 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CDC37
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215040.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q16543, Q61081, Q63692
Entrez	12539, 114562, 11140
Expiration Date	12 months from date of receipt.