

Product Datasheet

KAT7 peptide (orb218397)

Description	The peptide is used to block anti-KAT7 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	KAT7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215041.
Purity	> 85%
Source	Synthetic
Uniprot ID	O95251, Q810T5, Q5SVQ0
Entrez	217127, 303470, 11143
Expiration Date	12 months from date of receipt.