

Product Datasheet

CHEK2 peptide (orb218399)

Description	The peptide is used to block anti-CHEK2 antibody reactivity.
Reactivity	Human, Porcine
Tested Applications	Blocking
Target	CHEK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215043.
Purity	> 85%
Source	Synthetic
Uniprot ID	096017
Entrez	11200
Expiration Date	12 months from date of receipt.