

Product Datasheet

KLK8 peptide (orb218400)

Description	The peptide is used to block anti-Kallikrein 8 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	KLK8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215044.
Purity	> 85%
Source	Synthetic
Uniprot ID	O88780, Q61955, O60259
Entrez	11202, 259277
Expiration Date	12 months from date of receipt.