

Product Datasheet

MRPL3 peptide (orb218402)

Description	The peptide is used to block anti-MRPL3 antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	MRPL3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215046.
Purity	> 85%
Source	Synthetic
Uniprot ID	P09001
Entrez	11222
Expiration Date	12 months from date of receipt.