

Product Datasheet

RPL35 peptide (orb218403)

Description	The peptide is used to block anti-RPL35 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	RPL35
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215047.
Purity	> 85%
Source	Synthetic
Uniprot ID	P42766 , Q6ZWW7 , P17078
Entrez	11224 , 66489 , 296709
Expiration Date	12 months from date of receipt.