

Product Datasheet

GPR176 peptide (orb218406)

Description	The peptide is used to block anti-GPR176 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GPR176
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215050.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q80WT4, Q64017, Q14439
Entrez	117257, 11245, 381413
Expiration Date	12 months from date of receipt.