

Product Datasheet

SLC6A14 peptide (orb218408)

Description	The peptide is used to block anti-SLC6A14 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	SLC6A14
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215052.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UN76
Entrez	11254
Expiration Date	12 months from date of receipt.