

Product Datasheet

RRAS2 peptide (orb218412)

Description	The peptide is used to block anti-RRAS2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RRAS2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215056.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62070, P62071
Entrez	66922, 22800
Expiration Date	12 months from date of receipt.