

Product Datasheet

SV2C peptide (orb218421)

Description	The peptide is used to block anti-SV2C antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SV2C
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215065.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q496J9 , Q9Z2I6 , Q69ZS6
Entrez	75209 , 29643 , 22987
Expiration Date	12 months from date of receipt.