

## Product Datasheet

### HIC2 peptide (orb218426)

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-HIC2 antibody reactivity. |
| <b>Reactivity</b>          | Human, Mouse, Rat                                           |
| <b>Tested Applications</b> | Blocking                                                    |
| <b>Target</b>              | HIC2                                                        |
| <b>Form/Appearance</b>     | Lyophilized powder                                          |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                |
| <b>Note</b>                | For research use only                                       |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb215070.    |
| <b>Purity</b>              | > 85%                                                       |
| <b>Source</b>              | Synthetic                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q96JB3</a> , <a href="#">Q9JLZ6</a>             |
| <b>Entrez</b>              | <a href="#">58180</a> , <a href="#">23119</a>               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                             |