

Product Datasheet

PIP5K1C peptide (orb218435)

Description	The peptide is used to block anti-PIP5K1C antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PIP5K1C
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215079.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5I6B8 , O70161 , O60331
Entrez	314641 , 18717 , 23396
Expiration Date	12 months from date of receipt.