

Product Datasheet

MAPK8IP2 peptide (orb218445)

Description	The peptide is used to block anti-JIP2 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MAPK8IP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215089.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9ERE9, Q13387
Entrez	60597, 23542
Expiration Date	12 months from date of receipt.