

Product Datasheet

SMUG1 peptide (orb218448)

Description	The peptide is used to block anti-SMUG1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SMUG1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215092.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6P5C5, Q53HV7, Q811Q1
Entrez	71726, 315344, 23583
Expiration Date	12 months from date of receipt.