

Product Datasheet

SSBP2 peptide (orb218450)

Description	The peptide is used to block anti-SSBP2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SSBP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215094.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CYZ8, P81877
Entrez	66970, 23635
Expiration Date	12 months from date of receipt.