

Product Datasheet

SLC24A2 peptide (orb218456)

Description	The peptide is used to block anti-NCKX2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SLC24A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215100.
Purity	> 85%
Source	Synthetic
Uniprot ID	O54701, Q9UI40
Entrez	84550, 25769
Expiration Date	12 months from date of receipt.