

Product Datasheet

ATXN10 peptide (orb218458)

Description	The peptide is used to block anti-Ataxin 10 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ATXN10
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215102.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UBB4, P28658, Q9ER24
Entrez	25814, 170821, 54138
Expiration Date	12 months from date of receipt.