

Product Datasheet

TNFAIP8 peptide (orb218459)

Description	The peptide is used to block anti-SCC-S2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	TNFAIP8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215103.
Purity	> 85%
Source	Synthetic
Uniprot ID	095379, Q921Z5
Entrez	25816, 106869
Expiration Date	12 months from date of receipt.