

Product Datasheet

BRMS1 peptide (orb218463)

Description	The peptide is used to block anti-BRMS1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	BRMS1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215107.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5M7T3, Q9HCU9, Q99N20
Entrez	25855, 293668, 100911841, 107392
Expiration Date	12 months from date of receipt.