

Product Datasheet

RPL36 peptide (orb218465)

Description	The peptide is used to block anti-RPL36 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RPL36
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215109.
Purity	> 85%
Source	Synthetic
Uniprot ID	P47964, Q9Y3U8, P39032
Entrez	58927, 25873
Expiration Date	12 months from date of receipt.