

Product Datasheet

POLR1A peptide (orb218467)

Description	The peptide is used to block anti-POLR1A antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	POLR1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215111.
Purity	> 85%
Source	Synthetic
Uniprot ID	054889, 035134, 095602
Entrez	83581, 25885, 20019
Expiration Date	12 months from date of receipt.