

Product Datasheet

GORASP2 peptide (orb218469)

Description	GORASP2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GORASP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215113.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99JX3 , Q9R064 , Q9H8Y8
Entrez	26003 , 70231
Expiration Date	12 months from date of receipt.